## SORBACT® FOAM DRESSING

EN – Instructions for Use

SV – Bruksanvisning

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DA – Brugsanvisning

NO – Bruksanvisning

ET – Kasutusjuhised

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SK – Návod na použitie

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**AU/NZ** Bayport Brands Pty Ltd, 21 Bayport Court, Mornington, Vic 3931, Australia

**LB** Food and Drug Corporation Sal, Byblos Bank building, Rachid Karameh street, Verdun, Beirut P.O. Box 11-8182, Lebanon, Tel. +961-1-862000, [www.fdc.com.lb](http://www.fdc.com.lb)

**IR** NIBYKAN Golden Pharma Mannheim GmbH, Augustaanlage 32, 68165 Mannheim, Germany, Tel. +49 621 43179339

Payk Daru Tosseh, No. 03, Karimkhanzand Ave, 7th of Tir Sq, Tehran, Iran. 1585917114, Tel. +98 21 54549000

**Device Description** Sorbact® Foam Dressing is a sterile, bacteria and fungi binding foam dressing. It consists of a green Sorbact® wound contact layer combined with a polyurethane foam with a semi-permeable backing. Sorbact® Foam Dressing absorbs and retains exudate, thereby reducing the risk of maceration and enabling a moist wound environment. The semi-permeable film backing allows excess fluid to evaporate.

**Mode of Action** Sorbact® binds microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus species*, *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. These microorganisms are removed from the wound each time the dressing is changed.

**Intended Use** Sorbact® Foam Dressing is intended for use in management of clean, contaminated, colonized or infected wounds with low to moderate exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

**Contraindications** There are no known contraindications to the use of Sorbact® Foam Dressing. Attention is drawn to the warnings and precautions.

### **Warnings and Precautions**

If the treated condition deteriorates, fails to improve or if a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on patients with a known hypersensitivity to any of the dressing components.

The green wound contact layer should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing.

Should the dressing dry or adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not use in combination with fatty products, such as ointments, creams or solutions, as they decrease the binding of microorganisms.

Do not use with oxidizing agents such as hypochlorite solutions or hydrogen peroxide, as these can break down the absorbent component of the dressing.

Sorbact® Foam Dressing is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination.

Do not re-sterilize.

Do not use if peel pouch is damaged.

Prior to radiation therapy remove Sorbact® Foam Dressing. A new dressing can be applied following treatment.

Sorbact® Foam Dressing is MR safe.

**Instructions for Use** **1.** Prepare the wound according to local clinical practice. **2.** Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by at least 2 cm. **3.** Remove the dressing from the pouch using an aseptic technique. **4.** If cutting the dressing, use an aseptic technique. Discard any open and unused dressing. **5.** Apply the dressing. Ensure that the green side of the dressing comes into direct contact with the complete wound surface. **6.** Apply an appropriate dressing fixation. **7.** The dressing change frequency depends on exudate levels and overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.

**Storage and Disposal** Keep dry and away from sunlight. Disposal should be made according to local environmental procedures.

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