## SORBACT® FOAM GENTLE BORDER

EN – Instructions for Use

SV – Bruksanvisning

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DA – Brugsanvisning

NO – Bruksanvisning

ET – Kasutusjuhised

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SK – Návod na použitie

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**AU/NZ** Bayport Brands Pty Ltd, 21 Bayport Court, Mornington, Vic 3931, Australia

**LB** Food and Drug Corporation Sal, Byblos Bank building, Rachid Karameh street, Verdun, Beirut P.O. Box 11-8182, Lebanon, Tel. +961-1-862000, [www.fdc.com.lb](http://www.fdc.com.lb)

**IR** NIBYKAN Golden Pharma Mannheim GmbH, Augustaanlage 32, 68165 Mannheim, Germany, Tel. +49 621 43179339

Payk Daru Tosseh, No. 03, Karimkhanzand Ave, 7th of Tir Sq, Tehran, Iran. 1585917114, Tel. +98 21 54549000

**Device Description** Sorbact® Foam Gentle Border is a sterile, bacteria and fungi binding, self-adherent absorbent wound dressing. It consists of a green Sorbact® wound contact layer combined with polyurethane foam, soft silicone adhesive borders and a vapor permeable polyurethane backing film. Sorbact® Foam Gentle Border absorbs and retains exudate, thereby reducing the risk of maceration and enabling a moist wound environment. The water proof and breathable backing film provides protection against external contamination and allows excess fluid to evaporate. The skin friendly soft silicone adhesive borders enable atraumatic dressing changes. The dressing can be used in compression therapy.

**Mode of Action** Sorbact® binds hydrophobic microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus species*, *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. These microorganisms are removed from the wound each time the dressing is changed.

**Intended Use** Sorbact® Foam Gentle Border is intended for use in management of clean, contaminated, colonized or infected wounds with moderate exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

**Contraindications** There are no known contraindications to the use of Sorbact® Foam Gentle Border. Attention is drawn to the warnings and precautions.

### Warnings and Precautions

If the treated condition deteriorates, fails to improve or if a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on patients with a known hypersensitivity to any of the dressing components.

The green wound contact layer should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing.

Should the dressing dry or adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

Do not use with oxidizing agents such as hypochlorite solutions or hydrogen peroxide, as these can break down the absorbent component of the dressing.

Sorbact® Foam Gentle Border is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination.

Do not re-sterilize.

Do not use if peel pouch is damaged.

Prior to radiation therapy remove Sorbact® Foam Gentle Border. A new dressing can be applied following treatment.

Sorbact® Foam Gentle Border is MR safe.

**Instructions for Use** **1.** Prepare the wound according to local clinical practice. Ensure that the surrounding skin is clean and dry. **2.** Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by 1 cm to minimize the risk of maceration. **3.** Remove the dressing from the pouch using an aseptic technique. **4.** Remove the protective film from the wound contact side of the dressing and apply the dressing. Do not stretch. Ensure that the green wound contact layer comes into direct contact with the complete wound surface. **5.** Gently press the borders to the surrounding skin. If necessary the dressing can be adjusted and reapplied. **6.** The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.

**Storage and Disposal** Keep dry and away from sunlight. Disposal should be made according to local environmental procedures.

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