Popis produktu

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EN
<div>Device Description</div> <p>Sorbact® NPWT Wound Filler is a bacteria and fungi binding wound dressing. It consists of a green Sorbact® wound filler. Sorbact® NPWT Wound Filler conforms to the wound bed providing negative pressure distribution without in-growth of granulation tissue.</p>
<div>Mode of Action</div> <p>Sorbact® binds microorganisms, such as <i>Staphylococcus aureus</i> (including MRSA), <i>Streptococcus species</i>, <i>Escherichia coli</i>, <i>Pseudomonas aeruginosa</i> and <i>Candida albicans</i>, as shown <i>in vitro</i>. These microorganisms are removed from the wound each time the dressing is changed.</p>
<div>Intended Use</div> <p>Sorbact® NPWT Wound Filler is intended for use in conjunction with negative pressure wound therapy systems in management of clean, contaminated, colonized or infected wounds.</p>
<div>Contraindications</div> <p>There are no known contraindications to the use of Sorbact® NPWT Wound Filler. Attention is drawn to the warnings and precautions.</p>
<div>Warnings and Precautions</div> <p>If the treated condition deteriorates, fails to improve or a side effect is observed, consult a physician or an appropriate health care provider.</p> <p>Do not use on patients with a known hypersensitivity to any of the dressing components.</p> <p>The dressing should always be applied in direct contact with the wound to allow microorganisms to bind to the dressing. Should the dressing dry or adhere to the wound, moisten or soak the dressing with sterile saline to assist removal and avoid disruption of the healing wound.</p> <p>Do not use in combination with fatty products, such as vaseline gauze, as they may decrease the binding of microorganisms.</p> <p>Always ensure that Sorbact® NPWT Wound Filler does not come in direct contact with bone and tendons, exposed blood vessels, anastomotic sites, organs or nerves. Always make sure that no pieces of dressing material are left in the wound after removal.</p> <p>Sorbact® NPWT Wound Filler is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination.</p> <p>Do not re-sterilize.</p> <p>Do not use if peel pouch is damaged.</p> <p>Prior to radiation therapy remove Sorbact® NPWT Wound Filler. A new dressing can be applied following treatment.</p> <p>Sorbact® NPWT Wound Filler is MR safe.</p>
<div>Instructions for Use</div> <ol style="list-style-type: none">Prepare the wound according to local clinical practice. Select an appropriate dressing size for the wound. Ensure that the surrounding skin is clean and dry. Remove the dressing from the pouch using an aseptic technique. If cutting the dressing, discard any open and unused dressing. Moisten the dressing with sterile saline. Loosely fill the wound cavity and ensure that the dressing is in direct contact with the wound surface. If more than one piece of dressing is used, check that all pieces are in direct contact with each other to ensure an even distribution of negative pressure. Count the pieces of dressing used and record the total in the patient’s records. Ensure that the wound filler has no direct contact with intact skin. The pressure level should be set according to the physician’s instruction and be based on indication and wound condition. Proceed as directed for the NPWT system used. <p>Sorbact® NPWT Wound Filler should be changed within 48-72 hours.</p>
<div>Storage and Disposal</div> <p>Keep dry and away from sunlight.</p> <p>Disposal should be made according to local environmental procedures.</p>

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ABICO Medical AB
ABIGO Medical AB, Ekonomivägen 5, SE-436 33 Askim, Sweden