SORBACT® SUPERABSORBENT

EN – Instructions for Use

SV – Bruksanvisning

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DA – Brugsanvisning

NO – Bruksanvisning

ET – Kasutusjuhised

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SK – Návod na použitie

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AU/NZ Bayport Brands Pty Ltd, 21 Bayport Court, Mornington, Vic 3931, Australia

LB Food and Drug Corporation Sal, Byblos Bank building, Rachid Karameh street, Verdun, Beirut P.O. Box 11-8182, Lebanon, Tel. +961-1-862000, www.fdc.com.lb

IR NIBYKAN Golden Pharma Mannheim GmbH, Augustaanlage 32, 68165 Mannheim, Germany, Tel. +49 621 43179339

Payk Daru Tosseh, No. 03, Karimkhanzand Ave, 7th of Tir Sq, Tehran, Iran. 1585917114, Tel. +98 21 54549000

Device Description Sorbact® Superabsorbent is a sterile, bacteria and fungi binding, superabsorbent wound dressing. It consists of a green Sorbact® wound contact layer combined with a superabsorbent core and a white, water repellent backing, which prevents exudate strike through. Sorbact® Superabsorbent absorbs and retains exudate, thereby reducing the risk of maceration and enabling a moist wound environment. The high absorption capacity of Sorbact® Superabsorbent helps decrease dressing change frequency. The dressing can be used in compression therapy.

Mode of Action Sorbact® binds microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus species*, *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. These microorganisms are removed from the wound each time the dressing is changed.

Intended Use Sorbact® Superabsorbent is intended for use in management of clean, contaminated, colonized or infected wounds with high to excessive exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

Contraindications There are no known contraindications to the use of Sorbact® Superabsorbent. Attention is drawn to the warnings and precautions.

Warnings and Precautions

If the treated condition deteriorates, fails to improve or if a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on eyes, mucous membranes or in wound cavities. The dressing swells considerably after fluid absorption.

Do not use on arterial bleeds or heavily bleeding wounds.

Do not use on patients with a known hypersensitivity to any of the dressing components.

The green wound contact layer should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing.

Should the dressing dry or adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of microorganisms.

Sorbact® Superabsorbent is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination.

Do not re-sterilize.

Do not cut the dressing.

Do not use if peel pouch is damaged.

Prior to radiation therapy remove Sorbact® Superabsorbent. A new dressing can be applied following treatment.

Sorbact® Superabsorbent is MR safe.

Instructions for Use **1.** Prepare the wound according to local clinical practice. **2.** Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by at least 2 cm. **3.** Remove the dressing from the pouch using an aseptic technique. **4.** Apply the dressing. Ensure that the green side of the dressing comes into direct contact with the complete wound surface. **5.** Fixate with tape at dressing edges or cover with bandage. Ensure the dressing is allowed to expand under the fixation. **6.** The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. If saturated the dressing shall be changed. Should the clinical condition allow, the dressing can be left in place for up to 7 days.

Storage and Disposal Keep dry and away from sunlight. Store below 40°C/104°F. Disposal should be made according to local environmental procedures.

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