SORBACT® SURGICAL DRESSING

EN – Instructions for Use

SV – Bruksanvisning

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DA – Brugsanvisning

NO – Bruksanvisning

ET – Kasutusjuhised

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SK – Návod na použitie

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AU/NZ Bayport Brands Pty Ltd, 21 Bayport Court, Mornington, Vic 3931, Australia

LB Food and Drug Corporation Sal, Byblos Bank building, Rachid Karameh street, Verdun, Beirut P.O. Box 11-8182, Lebanon, Tel. +961-1-862000, www.fdc.com.lb

IR NIBYKAN Golden Pharma Mannheim GmbH, Augustaanlage 32, 68165 Mannheim, Germany, Tel. +49 621 43179339

Payk Daru Tosseh, No. 03, Karimkhanezand Ave, 7th of Tir Sq, Tehran, Iran. 1585917114, Tel. +98 21 54549000

Device Description Sorbact® Surgical Dressing is a sterile, bacteria and fungi binding wound dressing. It consists of a green Sorbact® wound contact layer combined with an absorbent wound pad and a transparent acrylic adhesive film. Sorbact® Surgical Dressing absorbs and retains exudate and enables a moist wound environment. The water proof backing film provides protection against external contamination and allows excess fluid to evaporate.

Mode of Action Sorbact® binds microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus species*, *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. These microorganisms are removed from the wound each time the dressing is changed.

Intended Use Sorbact® Surgical Dressing is intended for use in management of clean, contaminated, colonized or infected dry to low exuding wounds, such as surgical wounds and traumatic wounds, e.g. lacerations, cuts and abrasions.

Contraindications There are no known contraindications to the use of Sorbact® Surgical Dressing. Attention is drawn to the warnings and precautions.

Warnings and Precautions

If the treated condition deteriorates, fails to improve or a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on patients with a known hypersensitivity to any of the dressing components.

The green wound contact layer should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing.

Should the dressing adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

In common with all adhesive products, cases of irritation to the surrounding skin may occur. It should also be noted that inappropriate use or too frequent dressing changes, particularly in patients with fragile skin, may result in superficial skin trauma.

Sorbact® Surgical Dressing is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination.

Do not re-sterilize.

Do not use if peel pouch is damaged.

Prior to radiation therapy remove Sorbact® Surgical Dressing. A new dressing can be applied following treatment.

Sorbact® Surgical Dressing is MR safe.

Instructions for Use **1.** Prepare the wound according to local clinical practice. Ensure that the surrounding skin is clean and dry. **2.** Select an appropriate dressing size for the wound. **3.** Remove the dressing from the pouch using an aseptic technique. **4.** Remove the protective film from the wound contact side of the dressing and apply the dressing. Ensure that the green wound contact layer comes into direct contact with the complete wound surface. **5.** Press the borders to the surrounding skin and remove the stabilizing film on top. Start at the blue marker. **6.** The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.

Storage and Disposal Keep dry and away from sunlight. Disposal should be made according to local environmental procedures.

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