

Dialkylcarbamoyl chloride-impregnated dressing for the prevention of surgical site infection in women undergoing cesarean section: a pilot study

Authors:

Stanirowski PI, Kociszewska A, Cendrowski K and Sawicki W

Paper:

Arch Med Sci. 2016;12(5):1036-1042.

Material and methods:

This was a single-blinded randomized, controlled pilot study that evaluated the presence of superficial and deep SSIs during the first 14 days after a cesarean section. Enrolled patients were randomly allocated to receive treatment with either a DACC impregnated dressing (Sorbact® Surgical Dressing) or a standard surgical dressing. The primary study outcome was development of superficial or deep SSI within the first 14 days after a cesarean section. The surgical team was blinded to the dressing type until the time of skin closure.

Introduction:

Incisional surgical site infections (SSIs) occur in approximately 1.8–9.2% of patients subject to a cesarean section and contribute to prolonged hospitalization time and increased treatment costs.

Objective:

The objective of this study was to compare the new DACC dressing with a standard surgical dressing in the prevention of SSI in women undergoing cesarean section.

Main activities/study assessments

At 48 h	Wound evaluation 1	<ul style="list-style-type: none">• Dressings removed and a first clinical wound evaluation performed (according to CDC criteria)• The presence of surgical wound dehiscence (separation of the skin, subcutaneous tissue and/or fascia resulting from an incisional SSI) evaluated• For patients with signs of wound infection the following was noted: the time the first symptom occurred, need for systemic antibiotic treatment and need for a surgical intervention• From patients with clinical symptoms of SSI, wound swabs for culture were collected
Day 3	Patient discharged	<ul style="list-style-type: none">• Patient informed about the need to report at hospital in case of symptoms suggesting a possible infection (fever, pus from the surgical site, redness, edema, warmth, pain or tenderness of the surgical site area)
Day 7	Revisit to hospital Wound evaluation 2	<ul style="list-style-type: none">• Skin sutures were removed• A second clinical wound evaluation (according to CDC criteria)• The presence of surgical wound dehiscence evaluated (see details in wound evaluation 1)• For patients with signs of wound infection (see details in wound evaluation 1)• From patients with clinical symptoms of SSI, wound swabs for culture were collected
Day 14	Revisit to hospital Wound evaluation 3	<ul style="list-style-type: none">• A third clinical wound evaluation (according to CDC criteria)• The presence of surgical wound dehiscence evaluation (see details in wound evaluation 1)• For patients with signs of wound infection (see details in wound evaluation 1)• From patients with clinical symptoms of SSI, wound swabs for culture were collected

Results:

- 162 women undergoing elective or emergency CS were enrolled in the study and randomly allocated to receive treatment; 81 patients per group.
- No significant differences between the two study groups regarding basic demographic characteristics.
- 20 patients did not report for the scheduled follow-up visits after 7 and 14 days and were excluded from the analysis, leaving a total of 142 patients with follow-up data; 71 patients per group.
- The rate of surgical site infections was 2.8% in the Sorbact® group and 9.8% in the control group ($p = 0.08$).
- The rate of wound infection requiring systemic antibiotic treatment was 7.0% in the control group whereas systemic antibiotic treatment was not necessary in the study group ($p = 0.03$).
- The pre-pregnancy body mass index (BMI) was found to be a significant predictor for SSI after cesarean section.

Conclusion:

This is the first published study of Sorbact® in wound infection prophylaxis after cesarean section and the results indicate a decreasing tendency (NS/0.08) of the SSI rate after use of Sorbact® dressings.

Patients with SSIs who received a standard surgical dressing required systemic antibiotic therapy significantly more frequently ($p = 0.03$).

Summary:

The study did not show a statistical significant difference related to its primary objective.

A clear trend ($p = 0.08$) for a decreasing SSI tendency was however observed. This observation was supported by a significantly less need for systemic antibiotics in the Sorbact® group.